

Streptococcus canis MLST scheme

Primers used for amplification and sequencing

Primer name	Sequence (5' - 3')	Sequence length (bp)	Reference (PMID)
gki-fwd gki-rev	GCAGATTTTATTGGTATCGGTATGG GTGAACGTAAGAATTCTCCTGCTG	498	19386831
gtr-fwd gtr-rev	GGAATTGATTTAAACATTATGTCAGGAG CACAATAACGCCGCCATCCATA	450	19386831
murI-fwd murI-rev	TTATGGTCCAAGACCTGCTGAGC TTTCAGGACTTGCCGTTGTGTAATA	438	19386831
mutS-fwd mutS-rev	AGGTCAGATGTTAGAGGCTAGG CCTAGTTCATCAAATAGAATAAGGGAA	405	19386831
recP-fwd recP-rev	TGTCCGCACCCTATCAATGGAT CATCTTTCACAAGGATATGTTGCC	459	19386831
Xptgc-up Xptgc-dn	TTACTTGAAGAACGCATCTTA ATGAGGTCACTTCAATGCCC	450	20668530
yqiZ-fwd yqiZ-rev	CAGATGCTTTTAACAATTACCACATGG CCCATTACATTCACGATTTCTGG	434	19386831